## **TECHNICAL DATA**



# Nature Fresh® Odor Controller

## SUGGESTED APPLICATIONS

- **Patient Rooms**
- Restrooms
- **Hotel Rooms**
- Storage Containers/Luggage
- Offices
- **Vehicles**
- **Student Dorms**
- **Dressing Rooms**
- **Boats**
- **Utility Closets**



## **Restoration Tip!**

Nature Fresh odor neutralizers provide an effective solution for eliminating smoke odors from books, documents, clothing, bedding, and more! Simply put the contaminated articles in a box or plastic bag, place one of the barrels in the container, and seal it. Depending on the severity of damage and the number of items in the box/bag, treatment may require more than one barrel, may take from several hours to a few days, or may remain packed with goods during extended storage.

#### **PURPOSE**

Nature Fresh barrels are easy-to-use, portable, and passive systems. The 3" barrel treats 1,000 ft<sup>3</sup> while the 5" treats 2,000 ft<sup>3</sup>. Electricity and batteries are not required—Nature Fresh barrels rely on gentle air movement passing over the neutralizing element to constantly refresh the air. It uses Vaportek's patented spill-proof essential oil membranes to release an odor-neutralizing dry vapor into the air. Based on proven technology originally developed for use in medical facilities, Nature Fresh sets the standard for treating an area's atmosphere, structure, and contents simultaneously. The barrels eliminate the potential for dangerous spills and wet spots on floors, and contain no petroleum distillates, inert ingredients, or ozone-depleting chemicals.

## **DETAILS & USAGE**

Compact, durable, and lightweight housing. Can be used freestanding or as a wall mount. Lasts up to 60 days. To adjust the uniquely designed air vents on the top and bottom of the barrel, twist to open or close for precise vapor output. Nature Fresh barrels are sold by the dozen. For best results, store Nature Fresh barrels in a cool, dry location and keep away from extreme heat. Dry vapor output will be affected by humidity, airflow, and temperature. Ideal dry vapor release occurs at 50-80°F. Vapor release significantly increases in excess of 90°F. The passive system does not create moisture or ozone, and will not harm fabric, plants, rubber, adhesives, electrical components, or remove pigment from leather. Vaportek products use essential oils. If an allergic reaction occurs, please discontinue use immediately. Approximate Dimensions: 2"W and either 3.5"H or 5.5"H. Weight: less than 1 lb.

## **Reorder Information:**

- 3" Barrel (#90-8000) Neutral
- 3" Barrel (#90-8000-35) Winter
- 5" Barrel (#90-9000) Neutral

NOTE: Because we are committed to improving our products, Vaportek reserves the right to change product specifications and/or availability without notice or obligation.



P: 800.237.6367 | F: 262.246.5065

www.vaportek.com

www.facebook.com/vaportek in www.linkedin.com/vaportek Sussex, WI USA